



[→ Product Website](#)

Software for medical professionals

The RadiNET Pro software from EIZO enables the centralised and hence efficient quality management of all hospital monitors connected within the network. The browser-based application provides easy access to the information the user needs. The administrator can conduct or schedule quality control operations from anywhere in the hospital easily and promptly. Inventory information, e.g., the model name and usage time of monitors, can be managed in RadiNET Pro. In addition, information such as the connected computers and graphics boards, the medical facility and the installation location can also be recorded. Along with that, detailed information on each monitor can be consolidated into a report on request.

- ✓ Centralised quality management for monitors in the hospital network
- ✓ Easy-to-use browser-based applications in the facility's network
- ✓ Save time with remote calibration with the aid of monitor's own sensors
- ✓ Unified asset management for monitors, graphics boards and computers
- ✓ Manages five to 8,000 monitors; perfect addition to the RadiCS quality control software
- ✓ Reporting function with information on connected monitors, their usage and results

Technical Data

REQUIREMENTS (ADMINISTRATOR PC)

SUPPORTED WEB BROWSER	Microsoft Edge 79 or newer / Google Chrome™
MIN. RESOLUTION	1024x768

SERVER-REQUIREMENTS

Operating systems	Windows Server 2022 Standard, Windows Server 2019 Standard, Windows Server 2016 Standard, Windows Server 2012 R2 Standard, Windows 11 Pro / Enterprise, Windows 10 Pro / Enterprise (64-bit)
Databases	SQL Server 2022 Standard / Express Edition, SQL Server 2019 Standard/Express Edition, SQL Server 2016 Standard / Express Edition SP2
HARD DRIVE CAPACITY	150 GB (min.)
RAM	8 GB (min.)

GENERAL/OPERATION

MAX. NUMBER OF MANAGED PCs	1000 PCs
MAX. NUMBER OF MANAGED MONITORS	8000 monitors
SUPPORTED LANGUAGES	German, English, French, Japanese, Chinese

Find your EIZO contact:
EIZO AG - Switzerland
Moosacherstrasse 6, Au
8820 Wädenswil ZH
Phone +41 44 782 24 40
www.eizo.ch

All product names are trademarks or registered trademarks of EIZO Corporation in Japan and other countries or their respective companies. Copyright © 2024 EIZO Europe GmbH, Belgrader Str. 2, 41069 Mönchengladbach, Germany. All rights, errors and modifications reserved. Latest update: 12.05.2024